



SPARK^{Mi}TM

Redefine the
Measuring Experience



SPARKMi™

Redefine the Measuring Experience



Spark Mi looks like an ordinary tabletop mirror, but it is actually a measurement tool. Your client looks in the mirror wearing their glasses, you simply click, and Spark Mi captures the image.



Just the one picture is all that's needed for computerized measuring of all the required parameters, with precision accuracy.



And sunglasses pose no problem. Spark Mi clearly sees the client's pupils even through dark lenses, making measuring as accurate and easy as with clear lenses.



You can say goodbye to gadgets, rulers, and customer discomfort. Measuring for new glasses is now quick, easy, comfortable and precise.

Spark Mi – All the measurements you need:

- Spark Mi **automatically** measures near and far PD
- Mono PD
- FH
- Frame Box (A, B, ED)
- DBL
- BVD, Face Form Angle & Pantoscopic tilt



ALL FROM ONE PICTURE!

Spark Mi – Market-leading advantages



MEASURES EVEN WITH DARK SUNGLASSES!



NO GADGETS, CLIPS OR RULERS REQUIRED



ALSO SERVES AS AN ORDINARY STYLISH MIRROR



CONNECTS TO ANY PC



EXPORTS MEASUREMENTS TO ORDERING SYSTEMS



INTUITIVE AND EASY TO USE

Spark Mi is a revolutionary new measuring device that gives your clients greater comfort and a hi-tech experience, while it saves you time and guarantees precision.



SHAMIR

Shamir UK Ltd.
Bridge House, Saxon Way, Bar Hill,
Cambridgeshire, CB23 8AP

Tel: 01954 785100
E-mail: info@shamirlens.co.uk
Web: www.shamirlens.co.uk